

1 CORD & PLACENTA COLLECTION



STEP 1

Milk blood from the cord back into the placenta. Cut cord close to the placenta (**minimum of 7 inches of cord**) into the **Saline Bottle**. Place **Saline Bottle** into **Saline Bottle Bag**.



STEP 2

Place placenta in the sterile bag labeled **Bag 1**. Remove as much air as possible and seal shut.



STEP 3

Place **Bag 1** into **Bag 2**. Remove as much air as possible and seal shut.



STEP 4


Place **Bag 2** into **Bag 3**. Remove as much air as possible and seal shut. Return to Collection Kit. (If there is not a QR code label, record name, date and time on Bag 3.)

2 MATERNAL BLOOD COLLECTION



STEP 1

After delivery, use the vacutainer blood tubes (provided in the Maternal Collection Kit) to collect maternal blood. Filling 1 red tube and 2 purple tubes to capacity.



STEP 2

Ensure the provided QR code labels are attached to each tube and return the tubes to the absorbent sheet.



STEP 3


Place blood tubes in the biohazard bag provided. Remove as much air as possible and seal shut.



STEP 4

See other side for NanoCool Activation and Collection Kit Packing Instructions.

3 NANOCOOL ACTIVATION



STEP 1

Remove the NanoCool cooling engine lid from the box and place on a hard, flat surface with the silver-side facing down.


STEP 2

Using only a moderate amount of force, press straight down, with your thumb, on the white button. Using excessive force may cause the NanoCool to fail.

STEP 3

Depending on the ambient temperature, the NanoCool Logo will then turn blue, indicating that the cooling engine has been activated. This may take up to 3 minutes.

4 COLLECTION KIT PACKING SEQUENCE



STEP 1

Activate Instant Cold Pack (see instructions on pack) and place flat at the bottom of the kit.

STEP 2

Place Bag 3 (holding the triple-bagged placenta) flat on top of the cold pack. If the maternal blood draw is delayed, place Bag 3 in fridge or on wet ice (do not freeze) until blood draw is complete.

STEP 3

Place the saline bottle bag and biohazard blood draw bag with vials upright to the side of Bag 3.

STEP 4

After activation, return the NanoCool lid to original position (**silver-side up**) and press firmly to ensure snug fit on the top of the kit.

STEP 5

After provider paperwork has been completed and signed by provider, place paperwork in the clear pouch on kit lid and seal pouch before closing kit.

Once the kit is packed, notify CellSure ASAP at (866) 646-6171



Scan QR code to watch a kit packing video